



FEATURES

- PicoFrit® Integrated Technology
- High-performance separation at nanoliter flow rates
- Reproducible results: tip-to-tip, injection-to-injection, column-to-column, lab-to-lab, study-to-study
- Fused-Silica or nanoViper™ compatible connections
- Easy to use, reproducible robust universal connectivity
- Portable - compatible with multiple platforms

INTEGRATED PICOFRIT® TECHNOLOGY

- PicoTip Powered™ - built on SilicaTip™ technology delivering the standard for high-performance and reproducibility
- Zero-Dead Volume – featuring PicoFrit® Integrated Technology, the revolutionary emitter and column in one
- Digital Column Packing - custom and digital column packing featuring PicoClear™
- Integrated Packaging - simple, scalable, easy-to-use solution enabling standardized, portable and unified workflows

SPECIFICATIONS

Nanobore Format:	75 µm x 15 µm Tip ID
Bed Length:	25 cm
Chemistry:	Reprosil-PUR 1.9µm C18 (120 Angstrom)
Temperature:	50° C
Connectivity:	Universal
Pressure Limit:	650 bar (9,400 psi)