



Standard Products

Product Category	Part Number	Product Description
PicoFrit® Pioneer™	PF360-75-10-N-5	PicoFrit Pioneer: 360 µm OD, 75 µm ID, 10 µm tip ID, uncoated, box of 5, 50 cm length, self-pack column
	PF360-75-15-N-5	PicoFrit Pioneer: 360 µm OD, 75 µm ID, 15 µm tip ID, uncoated, box of 5, 50 cm length, self-pack column
	PF360-50-10-N-5	PicoFrit Pioneer: 360 µm OD, 50 µm ID, 10 µm tip ID, uncoated, box of 5, 50 cm length, self-pack column
SilicaTip™	FS360-20-10-N-20	SilicaTip: 360 µm OD, 20 µm ID, 10 µm tip ID, uncoated, box of 20, 5 cm length
	FS360-20-10-N-20-7CT	Precut SilicaTip: 360 µm OD, 20 µm ID, 10 µm tip ID, uncoated, box of 20, 7 cm length
	FS360-20-10-N-20-105CT	Precut SilicaTip: 360 µm OD, 20 µm ID, 10 µm tip ID, uncoated, box of 20, 10.5 cm length
	FS360-20-10-N-20-C12	SilicaTip: 360 µm OD, 20 µm ID, 10 µm tip ID, uncoated, box of 20, 12 cm length
	FS360-20-10-N-5	SilicaTip: 360 µm OD, 20 µm ID, 10 µm tip ID, uncoated, box of 5, 5 cm length
	FS360-20-10-N-5-7CT	Precut SilicaTip: 360 µm OD, 20 µm ID, 10 µm tip ID, uncoated, box of 5, 7 cm length
	FS360-20-10-N-5-105CT	Precut SilicaTip: 360 µm OD, 20 µm ID, 10 µm tip ID, uncoated, box of 5, 10.5 cm length
	FS360-20-10-N-5-C12	SilicaTip: 360 µm OD, 20 µm ID, 10 µm tip ID, uncoated, box of 5, 12 cm length
GlassTip™	BG12-94-4-CE-20	GlassTip: 1.2 mm OD, 0.94 mm ID, 4 µm tip ID, coated, box of 20
Accessories	Scribe	Diamond wedge scriber

Please inquire about the availability of custom sizes.