


# SOLUTIONS™

## SOLUTIONS FOR LIFE




UFS shown with optional equipment

### FEATURES

- Universal platform for nanoflow and microflow LC-MS operations
- Supports all PicoChip® and FlowChip™ consumables
- Effortless transition from nanoflow to microflow
- Integrated high-voltage, sheath gas and temperature control
- Tool-free installation for PicoChip and FlowChip consumables
- Push button operation to engage voltage
- Reproducible tip positioning
- Automatic magnetic alignment of chips
- Single connection to LC system
- Protected source environment
- Supports external trapping column (optional)
- USB Video microscope visualization (optional)
- Simple to install
- Easy to use—click and collect

### SPECIFICATIONS

|                       |   |
|-----------------------|---|
| MS Systems Supported: | Thermo Ion-Max™ and Ion-Max™ Next Gen equipped MS systems |
| ESI Voltage:          | MS control and delivery (0.3 to 3.5 kV typical)           |
| Sheath Gas:           | MS control and delivery                                   |
| Column Flow Rates:    | 200 nL/min to 10 µL/mi                                    |
| Column Temp:          | 35 to 65C   |
| XYZ Range:            | 6 mm each axis  |
| Video Microscope:     | PC host via USB (optional)                                |



USB Video Microscope (optional)



[newobjective.com](http://newobjective.com)  
[sales@newobjective.com](mailto:sales@newobjective.com)  
1 781 933 9560